COOK

MEDICAL

CE
0088

EN
6

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
Fig. 7


Fig. 8

NATHANSON TRANSCYSTIC BILE DUCT STONE EXPLORATION PACK

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Nathanson Transcystic Bile Duct Stone Exploration Pack is packaged preloaded and sterilized. The preloaded assembly includes: **(Fig. 1)**

- 5.5 French, 70 cm long radiopaque nylon catheter with a soft tip and five visual markers spaced 2 cm apart, starting at the tip of the catheter.
- 75 cm long, four-wire Flatwire Stone Extractor and plastic pin vise.
- "Y" Tuohy-Borst Adapter

This set is designed for use through a laparoscopic cholangiography instrument clamp (e.g., Reddick-Olsen Clamp) with a minimum 6 French working channel.

INTENDED USE

The Nathanson Transcystic Bile Duct Stone Exploration Pack is intended for extraction of bile duct stones via a transcystic approach.

CONTRAINDICATIONS

- Large bile duct stones which cannot be extracted through the cystic duct, as identified by evaluation of cystic duct diameter and stone diameter at the time of surgery

WARNINGS

- Rigidity of the catheter assembly may cause cystic duct and bile duct perforation during biliary cannulation, especially when acute inflammation is present or deployment past an impacted stone is encountered.
- Wire deployment within the ampulla of Vater and duodenum may entrap mucosae anchoring the catheter assembly *in situ*, potentially injuring mucosa and triggering pancreatitis.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in the anatomy of the biliary tree and stone extraction techniques.
- It is important that the surgeon and operating room personnel become completely familiar with the operation of this device prior to use.
- Adequate lead shielding for all operating room staff is required during fluoroscopic X-ray screening.
- Patient lead shielding is necessary if any suspicion of early pregnancy exists.
- Kinking of the catheter's soft tip may obstruct advancement of the Flatwire Stone Extractor.
- Stone extraction via the cystic duct may lead to stone fragmentation.
- Withdrawal of the Flatwire Stone Extractor back into the catheter, in preparation for biliary re-cannulation, may lead to entanglement in the laparoscopic cholangiography instrument clamp's tip mechanism, damaging the wires.
- Avoid opening the wire basket through or beyond the papilla. **(Fig 6)**

POTENTIAL ADVERSE EVENTS

- Fragmentation of stones in the bile duct
- Ampullary edema
- Bile duct perforation
- Pancreatitis
- Retained bile duct stones
- Retained stones in the peritoneal cavity
- Splitting of cystic duct

INSTRUCTIONS FOR USE

Preparation for Use

1. Remove contents from the package in a sterile manner and place on a sterile field.
2. Examine the contents to ensure that all components are intact and assembled correctly.
3. Withdraw the Flatwire basket back into the catheter approximately 5 cm by grasping the pin vise and stiff end of the exposed wire and gently pulling backwards. Note the position of the pin vise in relation to the valve. Maintain this position to ensure that the tip of the basket is not exposed at the distal end of catheter when inserted into the patient.
4. Gently tighten the Tuohy-Borst adapter to ensure that flush does not escape and to stabilize the wire basket for insertion. Do not over-tighten on the wire.
5. Flush the catheter system with sterile saline through the side-arm adapter to remove any air bubbles.
6. Load the assembly into a pre-positioned, open-tipped laparoscopic cholangiography clamp, ready for use.

7. Loosen the Tuohy-Borst adapter only after the catheter has passed the stone, and the wire and basket need to be advanced for stone capture.

Access and Cholangiography

1. Initiate standard laparoscopic exposure with umbilical and three subcostal ports.
2. During dissection of the cystic duct, note the duct size and dissect down closer to the common hepatic duct to avoid spiral valves.
3. Using an appropriately sized open-ended cholangiography catheter under fluoroscopic imaging, note the cystic duct course into the CBD (direct, into the right hepatic duct, or very low down near the ampulla) and the CBD diameter, coupled with the size, shape and number of bile duct stones. **NOTE:** Identification of stones in the biliary tree proximal to the cystic duct entry may be reason to choose either choledochotomy or ERCP clearance.
4. If there is uncertainty in interpreting images, it may be helpful to wash contrast out with saline and re-run the cholangiogram. Imaging during the washout phase can be helpful, as during the contrast run, to demonstrate details of the stones.

Transcystic Stone Extraction

1. Once bile duct stones have been confirmed and are considered suitable for transcystic extraction, exchange the *in situ* cholangiography catheter for the pre-prepared Nathanson 5.5 French catheter with pre-loaded basket, loaded into an appropriate laparoscopic cholangiography instrument clamp. Please refer to the "Preparation Prior To Use" section of these Instructions for Use for information on preparing the catheter. **NOTE:** Using a larger-sized catheter for the diagnostic cholangiogram may save time for patients who are at high risk for bile duct stones.
2. Under fluoroscopic monitoring, advance the catheter tip just beyond the first stone, then extend the Flatwire basket down the catheter until its tip emerges. (**Fig. 3**)
3. Hold the basket's position steady and withdraw the catheter and Tuohy-Borst adapter until the adapter abuts the basket handle. This will allow the basket to deploy around the stone. (**Fig. 4**)
4. Engage the wires around the stone using a gentle vibrating motion. Upon capturing the stone, slowly withdraw the basket and catheter as a unit from the cystic duct. (**Fig. 5**)
5. Place the stone in a suitable area of the abdominal cavity for its later removal, along with the gallbladder. **NOTE:** Calculi may be fragmented upon removal. Further visual inspection with a small-caliber, flexible choledochoscope may be helpful to ensure complete clearance of the stone.
6. In the presence of multiple stones, repeat the process in sequence, working down toward the ampulla. **CAUTION:** Avoid deploying the basket in the ampulla of the distal common duct, as its mucosa may catch in the basket wires, potentially inducing acute pancreatitis or resulting in long-term ampullary scarring. Early resistance upon basket withdrawal and fluoroscopic inversion of the ampulla and the distal bile duct are indicators of entanglement. (**Fig. 7**) Should this occur, discontinue catheter withdrawal and disentangle the ampulla mucosa from the basket wires by pushing the basket and the tip of the catheter into the duodenum. (**Fig. 8**) Withdraw the basket into the catheter, and withdraw the catheter into the bile duct.
7. For an impacted stone, pass a 5 French balloon catheter past the stone, inflate the balloon and, by traction, displace the stone. **NOTE:** Once the stone is sufficiently dislodged and contrast can flow past it, retrieve the stone with the Flatwire basket, as ongoing balloon extraction may result in stone displacement into the proximal bile ducts.

HOW SUPPLIED

Supplied Sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature.

Taylor CJ, Kong J, Ghusn M, White S, Crampton N, Layani L. Laparoscopic bile duct exploration: Results of 160 consecutive cases with 2-year follow-up. *ANZ J. Surg.* 2007; 77: 440-445

Nathanson LK, O'Rourke NA, Martin IJ, Fielding GA, et al. Postoperative ERCP versus laparoscopic choledochotomy for clearance of selected bile duct calculi: a randomized trial. *Annals of Surgery* 2005; 242: 188-192

Sgourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extration and laparoscopic cholecystectomy for choledocholithiasis. A prospective randomized study. *Minerva Chirurgica* 2002; 57: 467-474

Urbach DR, Khanjanchee YS, Jobe BA et al. Cost-effective management of common bile duct stones. *Surg Endosc* 2001; 15: 4-13

Liu TH, Consorti ET, Kawashima A, Tam EP, et al Patient Evaluation with selective use of Magnetic Resonance Cholangiography and Endoscopic Retrograde Cholangiopancreatography Before Laparoscopic Cholecystectomy. Ann Surg 2001; 234: 33-40

Traverso LW. A cost analysis of the treatment of common bile duct stones discovered during cholecystectomy. Semin Laparosc Surg 2000; 7: 302-307

Cuschieri A, Lezoche E, Mornino M, Croce E, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management of patients with gallstone disease and ductal calculi. Surgical Endoscopy 1999; 13: 952-957

Rhodes M, Sussman L, Cohen L, Lewis MP. Randomized trial of laparoscopic exploration of common bile duct versus postoperative endoscopic retrograde cholangiography for common bile duct stones. Lancet 1998; 351: 159-161

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Taylor CJ, Kong J, Ghosn M, White S, Crampton N, Layani L. Laparoscopic bile duct exploration: Results of 160 consecutive cases with 2-year follow-up. *ANZ J. Surg.* 2007; 77: 440-445

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Liu TH, Consorti ET, Kawashima A, Tam EP, et al Patient Evaluation with selective use of Magnetic Resonance Cholangiography and Endoscopic Retrograde Cholangiopancreatography Before Laparoscopic Cholecystectomy. *Ann Surg* 2001; 234: 33-40

Traverso LW. A cost analysis of the treatment of common bile duct stones discovered during cholecystectomy. *Semin Laparosc Surg* 2000; 7: 302-307

Cuschieri A, Lezoche E, Mornino M, Croce E, et al. E.A.E.S. multicenter prospective randomized trial comparing two-stage vs single-stage management of patients with gallstone disease and ductal calculi. *Surgical Endoscopy* 1999; 13: 952-957

Rhodes M, Sussman L, Cohen L, Lewis MP. Randomized trial of laparoscopic exploration of common bile duct versus postoperative endoscopic retrograde cholangiography for common bile duct stones. *Lancet* 1998; 351: 159-161

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Taylor CJ, Kong J, Ghusn M, White S, Crampton N, Layani L. Laparoscopic bile duct exploration: Results of 160 consecutive cases with 2-year follow-up. ANZ J. Surg. 2007; 77: 440-445

Nathanson LK, O'Rourke NA, Martin IJ, Fielding GA, et al. Postoperative ERCP versus laparoscopic choledochotomy for clearance of selected bile duct calculi: a randomized trial. Annals of Surgery 2005; 242: 188-192

Sgourakis G, Karaliotas K. Laparoscopic common bile duct exploration and cholecystectomy versus endoscopic stone extration and laparoscopic cholecystectomy for choledocholithiasis. A prospective randomized study. Minerva Chirurgica 2002; 57: 467-474

Urbach DR, Khanjanchee YS, Jobe BA et al. Cost-effective management of common bile duct stones. Surg Endosc 2001; 15: 4-13

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MANUFACTURER
COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

EC REPRESENTATIVE
WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632
Bjaeverskov, DENMARK