Fanelli Laparoscopic Endobiliary Stent System

Instructions for Use

**EN
2-5**

Fanelli laparoskopisches endobiliäres Stentsystem

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**DE
6-9**

Sistema De Stent Endobiliar Laparoscópico De Fanelli

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**ES
10-13**

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**IT
18-21**

Fanelli Laparoscopic Endobiliary Stent System

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Eliminates the need for external drainage catheters and T-tubes typically required after common bile duct exploration. It is used to provide internal drainage in procedures where all stones have been cleared from the common bile duct. It may also be useful in patients with cholangitis, in whom surgical therapy is undertaken rather than ERCP. It is ideally suited for placement when cholangiography demonstrates common bile duct stones that are not removed at the time of cholecystectomy. The stent system is designed to securely maintain the stent within the delivery mechanism until the surgeon actively deploys the stent, allowing the surgeon to maintain control over the stent at all times until it is deployed. The delivery system features a sideport that invites contrast injection throughout the placement process, even while the wire guide is in place. Four integrated radiographic markers facilitate positioning of the stent, and make stent position clear to the surgeon throughout placement. These markers also visually confirm proper deployment of the stent, and inform the surgeon when the delivery mechanism has cleared the stent, and can be removed. The design also permits performance of a completion cholangiography to document stent placement and verify common bile duct drainage.

INTENDED USE

The Fanelli Laparoscopic Endobiliary Stent System is designed for internal biliary drainage. The product is intended for use by physicians trained and experienced in general surgical techniques.

CONTRAINDICATIONS

None Known

WARNINGS

Once the outer sheath has been retracted, it may not be possible to retrieve the stent using the delivery mechanism.

PRECAUTIONS

- It is important that the surgeon and operating room personnel become completely familiar with the operation of this device prior to first use.
- The stent should be advanced only over a wire guide, and should never be advanced through the common bile duct or into the duodenum across the ampulla, without a wire guide in place.
- Accurate placement of the stent requires fluoroscopic imaging, and it should not be used in the setting of static cholangiography.

INSTRUCTIONS FOR USE

The following illustration describes the components of the stent delivery system. Please refer to the accompanying illustration and study the position and relationship of the parts of this system. The inner carrier (a) supports the stent, and bears the radiographic markers that are important visual clues during stent placement. The stent (b) is the device that will be left across the ampulla, draining the common bile duct; it is a 7 French double flange stent. The stent backstop (c) is the backbone of the system, and serves the important function of maintaining the stent in proper position during deployment. The backstop can also be used to perform completion cholangiography if desired. The outer sheath (d) is the system component that prevents premature deployment of the stent. It secures the stent to the delivery mechanism, and must be actively withdrawn at the time of stent deployment, making accidental deployment unlikely.


Fig. 1

The purple (e) and white (f) Luer connectors represent the second step in active deployment of the stent. This color-coded assembly permits the surgeon to withdraw the inner carrier successfully, even in low light conditions, to complete deployment of the stent. The terminal two-way adapter (g) is designed to accept a wire guide smoothly, even when back loading the device onto a wire already in place across the ampulla. It also permits contrast injection throughout the procedure to maintain optimal opacification of the common bile duct and duodenum, ensuring accurate stent placement. The 150 cm .035 inch Tracer Hybrid wire guide (h) completes the system, serving as the guide to delivery of the stent assembly through the common bile duct and into the trans-ampullary position. The sealed introducer sheath permits placement of the 8.5 French stent assembly through a 10-12 mm laparoscopic trocar, or through a separate puncture site, depending on surgeon preference. The dilation balloon permits dilation of the cystic duct to allow placement of the stent via the transcystic route.

1. Remove the delivery mechanism and stent from the peel packs using sterile technique.
2. Familiarize yourself with the components of the system; inner carrier, stent, stent backstop, outer sheath, Luer connectors and two-way terminal adapter.
3. Be sure the junction of the purple (e) and white (f) Luer connectors is secure, but only finger tight.
4. At this point, the 150 cm .035 inch Tracer Hybrid wire guide (h) should be advanced through the wire port of the cholangiography catheter, and directed fluoroscopically through the catheter, distal common bile duct, and across the ampulla to coil gently within the duodenum.


Fig. 2

5. About 15-20 cm of the flexible wire guide should be coiled within the duodenum to prevent inadvertent removal of the wire guide, and the position of the wire should be confirmed fluoroscopically. (Fig. 2)

ENGLISH

6. The cholangiogram catheter and cholangiography clamp are then removed over the Tracer Hybrid wire guide, maintaining good wire position across the ampulla.
7. Introduction of the Fanelli Laparoscopic Endobiliary Stent System differs slightly depending on the method used to perform cholangiography.
8. If the cholangiogram catheter was introduced through:
 - a. A standard 10-12 mm epigastric trocar using a cholangiogram clamp, the sealed reducing sheath should be introduced over the wire to size the trocar appropriately for stent placement.
 - b. A standard 5 mm epigastric trocar using a cholangiogram clamp, no sheath or reducer is necessary and the stent delivery device will be introduced directly over the wire.
 - c. A separate puncture site using an angiocath, the cholangiogram catheter will be removed over the wire guide, followed by removal of the angiocath, so that the wire guide remains in place across the ampulla. The sealed reducing sheath, complete with its tapered tip dilator, is then introduced over the wire after enlarging the angiocath incision slightly. The dilator is then removed while the wire guide remains in place for introduction of the stent.
9. The stent may be placed via the cystic duct or directly through a choledochotomy.
10. Before introducing the stent system using the transcystic route, a visual assessment of the size of the cystic duct is necessary.
11. The outer diameter of the stent system is 8.5 French, and it is advisable to dilate the cystic duct using an over the wire biliary dilation balloon so the duct will smoothly accommodate introduction of the stent system.
12. The Fanelli Laparoscopic Endobiliary Stent System will then be introduced over the wire guide, under fluoroscopic control. Tightening of the Tuohy-Borst valve during injection of contrast will minimize leakage.


Fig. 3

13. As the stent system is introduced, the surgeon will detect four radiographic markers during fluoroscopic visualization, and the side access port can be used for contrast injection at any time additional visual enhancement is necessary. **(Fig. 3)**
14. The first radiographic marker (1) appears at the distal tip of the inner guide, matching the position of the distal tip of the stent itself.
15. The second marker (2) appears on the inner guide at a position matching the distal flange of the stent, and indicates the proper position of the stent across the ampulla. When this marker is seen in the duodenum, the surgeon is assured that the stent will be deployed across the ampulla, with the distal flange within the duodenum and the proximal flange within the common bile duct.

16. The third marker (3) denotes the position of the proximal flange, to demonstrate the position of both flanges to the surgeon before deployment. This marker should be kept above the ampulla to avoid deployment of the stent in a position that is too low in the bile duct.
17. The fourth marker (4) is the largest, and represents the distal tip of the stent backstop (c), against which deployment occurs.
18. After the stent assembly has been directed across the ampulla over the wire guide, markers 1 and 2 should be positioned within the duodenum, and markers 3 and 4 within the bile duct. It is deployed as follows:
19. The outer sheath (d), which securely attaches the stent to the delivery mechanism, is slowly retracted back toward the purple (e) and white (f) Luer hubs, until it contacts the white hub. This is done with the surgeon supporting the Luer hub in position without movement. Frequent fluoroscopic evaluation for proper positioning is advised. The surgeon should check positioning of the radiographic markers frequently, to be sure the stent location has not changed.
WARNING: Once the outer sheath has been retracted, it may not be possible to retrieve the stent using the delivery mechanism.
20. The wire guide (h) is then removed completely.
21. The surgeon then unscrews the white Luer hub (f) from the purple (e) Luer connector, holding the purple connector securely without permitting it to move, while drawing the white hub away from the assembly. This is performed using continuous fluoroscopy to maintain appropriate positioning of the radiographic markers, thereby maintaining positioning of the stent.
22. When markers 3, 2, and finally 1 have been drawn past marker 4, the surgeon has visual confirmation that the stent has been successfully deployed, and is completely separated from the delivery mechanism.
23. Following deployment, the surgeon may obtain a completion cholangiogram by injecting contrast through the purple (e) Luer adapter before removal of the stent backstop (c). The backstop should be withdrawn into the proximal common bile duct before cholangiography is performed to avoid inadvertent contact with the previously placed stent.
24. Once the entire assembly has been removed, the cystic duct stump is ligated according to surgeon preference, and cholecystectomy completed in the typical manner.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.