COOK
MEDICAL

CE 0123

EN
3

Zenapro™
Hybrid Hernia Repair Device

Instructions for Use

DA
5

Zenapro™
Hybridprodukt til herniereparation

Brugsanvisning

DE
7

Zenapro™
Hybrid Hernienimplantat

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ES
9

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
ZENAPRO™ HYBRID HERNIA REPAIR DEVICE

INTENDED USE

The Cook® Zenapro™ Hybrid Hernia Repair Device is intended to be implanted to reinforce soft tissues where weakness exists. Indications for use include the repair of a hernia or body wall defect requiring the addition of a reinforcing or bridging material to obtain the desired surgical result. The device is supplied sterile and is intended for one-time use.

Rx ONLY This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Surgeon should be trained in management of complications resulting from these procedures.

Read this document carefully. Failure to follow instructions may cause malfunction and/or cause patient injury.

This document is not a training manual on the repair of a hernia or body wall defect. The Cook® Zenapro™ Hybrid Hernia Repair Device should only be used by surgeons instructed in the surgical correction of this pathology and in the use of this device.

HYBRID HERNIA REPAIR DEVICE This symbol means the following: Hybrid

Hernia Repair Device

This device is intended for use by trained medical professionals.

CONTRAINDICATIONS

This device is derived from a porcine source and should not be used in patients with known sensitivity to porcine material. The Hybrid Hernia Repair Device is contraindicated for use in infected wounds.

WARNINGS

When this device is used in infants, children, pregnant women or women planning pregnancies, the surgeon should be aware that this product will not stretch significantly as the patient grows.

PRECAUTIONS

- Use of this device in repair of large hernias or body wall defects has not been clinically evaluated.
- **This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.**
- **Do not resterilize.** Discard all open and unused portions of the device.
- The device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that device is rehydrated prior to suturing, stapling, tacking or loading of the device laparoscopically.
- Ensure that all layers of the device are secured when suturing, stapling or tacking.
- Place device in maximum possible contact with healthy, well-vascularized tissue to encourage cell ingrowth and tissue remodeling.
- Suturing, stapling or tacking more than one device together may decrease device performance.
- No studies have been conducted to evaluate the reproductive impact of the clinical use of the device.
- **Extended rehydration or excessive handling at device edges where synthetic material is not present could lead to delamination of the device. To reduce this likelihood when handling, positioning, or tensioning the graft, grasp the graft in areas where both the biological and synthetic layers are present.**
- Care should be taken to avoid damage to the device when loading laparoscopically. It is recommended to load through a 10 mm or larger port.
- **Use the device as supplied without cutting or trimming. When device is placed extraperitoneally, with no exposure to abdominal organs, cutting and trimming is permitted.**

POTENTIAL COMPLICATIONS

Possible adverse reactions with the use of any prosthesis may include, but are not limited to:

- infection
- inflammation
- allergic reaction

- adhesion
- erosion
- fistula formation
- seroma formation
- hematoma
- recurrence of tissue defect

Complications such as delayed wound infection, hernia recurrence and the need for re-operation should be reasonably expected in patients who are critically ill.

As with any surgical procedure, there are always risks of other complications. These include but are not limited to:

- wound dehiscence
- wound complications
- device contraction
- fever
- pain
- bowel obstruction
- ileus
- revision/resurgery^{1,2}

If conditions of infection, inflammation, or allergic reaction cannot be resolved, consider removal of the device.

STORAGE

This device should be stored in a clean, dry location at controlled room temperature.

STERILIZATION

This device has been sterilized with ethylene oxide.

INSTRUCTIONS FOR USE

REQUIRED MATERIALS

- Sterile forceps
- Rehydration fluid: room temperature sterile saline or sterile lactated Ringer's solution.

NOTE: Always handle the device using aseptic technique, minimizing contact with latex gloves.

1. Using aseptic technique, remove the sterile inner package from the outer packaging and place it in the sterile field.
2. Open the sterile inner package carefully. Aseptically remove the device from the inner package.
3. Rehydrate the device in room temperature sterile saline or sterile lactated Ringer's solution until the desired handling characteristics are achieved. Rehydration time of greater than 1 minute is not required prior to suturing, stapling, tacking or laparoscopic loading of the device.
4. Prepare the site using standard surgical techniques.

DEVICE PLACEMENT

- Preperitoneal/retrorectus (Rives-Stoppa technique) or intraperitoneal placement (underlay) is recommended for open procedures.
- Intraperitoneal placement (underlay) is recommended for a laparoscopic approach.

5. **When placing the device, care must be taken to orient the device properly. The rough side should be placed against the abdominal wall.**

6. Fundamental surgical principles suggest recurrence can be minimized if the device overlaps the surrounding tissue by at least 4-5 cm in all directions.³

NOTE: If the chosen device is too small for the defect, excess tension may be placed on the suture line. This can result in recurrence of the original tissue defect or development of a defect in the adjacent tissues.

7. Using aseptic technique, transfer the device to the surgical site and suture, staple or tack into place, avoiding excess tension.

NOTE: Excessive handling may lead to partial delamination of a superficial layer. This can be managed by securing the detached portion to adjacent tissue with sutures or staples.

NOTE: The liberal use of transfascial sutures is recommended. Tacking devices alone may not provide adequate fixation to prevent recurrence.

NOTE: Use of surgical tacks may disrupt the smooth surface of the device and should be used sparingly when the device is in contact with viscera.

NOTE: Surgical experience indicates that suturing or stapling the device with close tissue approximation produces better outcomes. Leave less than or equal to 3 cm between sutures with bite depth of 1.5 cm to ensure the blue mesh is captured. Use permanent or long-term absorbable sutures or tacks according to surgeon preference.

NOTE: Interrupted sutures can provide additional security against recurrence of tissue defect in the event of suture failure.

8. Complete the standard surgical procedure.
9. Discard any unused portions according to institutional guidelines for medical waste.

10. Place closed suction drains for 2-6 weeks. Remove when output is less than 20 mL/24 hours for at least two (2) consecutive days or until drain is dry.

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West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.