EN
3

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T _ A U B _ R E V 1

ASCEND® AQ® URETERAL DILATION BALLOON CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque balloon catheter
- Inflation device (*optional*)

INTENDED USE

Used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Radiopaque markers indicate the proximal and distal ends of the balloon.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.
- Do not overinflate. Using excessive pressure to inflate the balloon on this device can cause the balloon to rupture.

PRECAUTIONS

- Balloon dilation catheters are intended for use by physicians trained and experienced in techniques for balloon-catheter dilation.
- Do not exceed the maximum rated burst pressure (listed on label) for this balloon device.
- Do not preinflate the balloon.
- Refer to product label or the inflation check valve on the balloon device for appropriate balloon volume
- The burst pressure data was analyzed using factors for a one-sided tolerance to determine with 95% confidence that 99.9% of these balloons would not burst at or below the calculated rated burst pressure.

POTENTIAL ADVERSE EVENTS

- Complications that might occur during this procedure could include overinflation of the balloon which could result in trauma to the surrounding tissue.

INSTRUCTIONS FOR USING ASCEND® AQ® URETERAL DILATION BALLOON CATHETER BALLOON PREPARATION

1. Upon removal from the package, inspect the balloon catheter to ensure no damage has occurred during shipment.
2. Remove the protective sleeve from the balloon and discard. Lock a 10 ml syringe filled with 5 ml of dilute contrast medium onto the balloon lumen and apply negative pressure for 20 to 30 seconds. Release negative pressure on syringe allowing medium to draw into balloon. Detach the syringe leaving medium in the balloon lumen.

BALLOON CATHETER INTRODUCTION AND INFLATION

1. Under fluoroscopic control, cystoscopically pass a wire guide (up to .038 inch in diameter) the desired distance into the ureter beyond the area of planned dilation.
2. Activate the AQ hydrophilic coating by immersing the balloon and catheter in a sterile physiological solution for approximately 30 seconds.

3. Carefully advance the balloon catheter over the previously placed wire guide under fluoroscopic guidance utilizing the radiopaque markers to ensure proper positioning. Advance catheter so the entire balloon extends beyond the end of the cystoscope before inflating the balloon.
4. Prepare inflation device in standard manner and purge all air from syringe and tubing. To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended.
5. Attach the prepared inflation device to the balloon lumen, open the stopcock on the inflation device and inflate the balloon.
6. Close the stopcock on the inflation device for at least 30 seconds to maintain pressure and allow adequate dilation. Time, rather than excessive pressure, is the key factor in transluminal dilation.

BALLOON DEFLATION AND WITHDRAWAL

1. Aspirate the balloon completely.
2. When the balloon is completely deflated, remove the device by withdrawing the catheter slowly. Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

If your device did not come with a Cook Sphere Inflation Device then the below instructions are not applicable, if your device did come with a Cook Sphere Inflation Device please refer to the instructions below for proper use.

INFLATION DEVICE

DEVICE DESCRIPTION

The Cook Sphere Inflation Device is a one-piece, plastic, disposable inflation device with a lock lever design that controls the piston, a manometer, and a connecting tube with a male rotating adapter. An optional 3-way stopcock may also be enclosed for use during preparation of the device. The manometer measures pressures ranging from vacuum to gauge capacity; the gauge is marked in 1 atm increments. The gauge also has an inner scale of comparable PSI measurements. The accuracy of the manometer has been determined to be within 1 atm over the range.

INTENDED USE

The inflation device is recommended for use while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Use only liquid inflation media. Do not inflate with air.
- Always follow the manufacturer's directions accompanying the balloon dilation catheter for instructions for use, maximum balloon inflation pressure, precautions, and warnings for that device.

PRECAUTIONS

- Before use, inspect the device to verify that no damage has occurred during shipping and handling.
- Before use, ensure the connector tubing is completely free of air.

INSTRUCTIONS FOR USE

PREPARATION

Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked.

Unlock the piston by pushing the lock lever left. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

1. Prepare a solution of contrast medium and normal saline in a small sterile bowl. Check balloon catheter and contrast medium instructions for specific contrast mixture recommendations.
2. Orient the tubing downward into the contrast medium.
3. Push the release lever left and aspirate enough solution to fill the syringe. (Attach stopcock, if applicable.)
4. Hold the device upright to purge the air from the syringe and connecting tube. Tap the syringe lightly, if necessary, to remove all the air bubbles and fill the connecting tube completely.
5. Inspect the syringe and tubing (and stopcock, if applicable) to ensure that the device has been completely purged of air bubbles.
6. Adjust the syringe volume to the desired amount. If more contrast is needed, submerge the syringe tip into the basin of solution and aspirate. (Close stopcock, if applicable.)

ATTACHING THE INFLATION DEVICE TO THE BALLOON DILATION CATHETER

1. Prepare and test the balloon dilation catheter according to the manufacturer's directions for use.
2. If a separate syringe was used to prepare the balloon catheter, remove it. When a stopcock is installed on the end of the inflation device connecting tube, it should be opened and purged with contrast media from the inflation device to eliminate air. Create a fluid-fluid connection between the balloon and the stopcock or connecting tube (male rotating adapter) of the inflation device by placing a drop of contrast solution from the syringe into each hub.
3. Hand-tighten the hubs securely.

OPERATING INFLATION DEVICE

1. Release the lock lever and allow the piston to move forward into neutral position (0 atm).
2. To inflate the balloon, engage the lock lever, turn the palm grip on the piston clockwise slowly until the desired inflation pressure is reached. The lock lever maintains the increasing pressure.
3. To gradually deflate the balloon, turn the palm grip on the piston counterclockwise slowly until the desired deflation pressure is reached.
4. To rapidly deflate the balloon, push the lock lever left, releasing the piston, and pull back. Slide the lock lever back to lock, if desired.
5. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations. For one time use only. Do not reuse. Do not re-sterilize. Read directions prior to use.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.