EN
4

Endoscopic Curved Needle Driver

Instruction for Use

CS
7

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
Figure E

ENDOSCOPIC CURVED NEEDLE DRIVER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Set Contains (*set components may vary*)

- Stainless steel shaft with retractable jaw at distal tip

INTENDED USE

Used for internal suturing under direct vision during intrabdominal endoscopic procedures. The needle driver design affords precise placement and control of standard curved suture needles, similar to open abdominal suturing.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

This product must be cleaned immediately after use. Blood, mucus and other debris, if allowed to dry, can render the device inoperable.

INSTRUCTIONS FOR USE

Endoscopic Curved Needle Driver

1. Advance the needle driver through a suture introducer.
2. Secure the proximal end of the suture in the instrument's jaws. **(Figure A)**
3. Withdraw the needle driver from the suture introducer, allowing the suture needle to extend from the distal tip of the introducer. **(Figure B)**
4. Disengage the suture from the needle driver jaws.
5. While stabilizing the suture inside the introducer, advance the needle driver through the introducer again.
6. Secure the proximal end of the suture needle (at needle suture junction) in the driver jaws, keeping the needle curve parallel with the needle driver and suture introducer. **(Figure C)**
7. Pull the suture taut (not tight) while withdrawing the loaded needle driver into the introducer until the suture needle is inside the distal end of the introducer.
8. Cut off the excess suture, leaving approximately one inch or less exposed from the proximal end of the suture introducer.
9. Advance the loaded needle driver and suture introducer through the trocar port as a unit, leaving the suture inside the introducer channel. With a grasper placed through another trocar port, readjust the position of the suture needle in the needle driver jaws. **(Figure D)**
10. With the needle positioned for suturing **(Figure E)**, advance the needle through the tissue to be sutured.
11. Following needle placement through the desired tissue, the entire length of the suture can be brought into the abdominal cavity and the suture can then be equalized in length in preparation for tying.
12. For needle removal, grasp the suture within the jaws of the needle driver 2 cm from the needle, cut the suture, and withdraw the needle into the introducer channel.

13. Withdraw the needle driver, suture needle, and suture introducer all at once through the trocar port, leaving the trocar port in place.
14. The suture is now ready for intra-abdominal tying.

Disassembly of the Needle Driver

1. Unscrew the end cap at the rear of the handle.

NOTE: Slight squeezing pressure on the handle loop will aid in the removal of the end cap.

2. Pull the handle loop down and rotate forward to remove.
3. Remove the inner slide mechanism out the rear of the handle.

Assembly of the Needle Driver

1. Insert the inner slide mechanism through the rear of the handle and push as far forward as possible.
2. Be sure the clamping wedge on the inner slide mechanism faces down before inserting.
3. Insert the front of the handle loop into the slot and rotate up into position.
4. Replace the end cap and screw tight.

NOTE: Slight squeezing pressure on the handle loop will aid in tightening the end cap.

Cleaning/Disinfections/Sterilization Instructions for the Endoscopic Curved Needle Driver

All cleaning/disinfection/sterilization methods listed herein have been validated.

Cleaning Instructions

1. Disassemble the needle driver per the instructions listed within the IFU.
 2. Remove all residual organic matter, blood and irrigation solution to allow the sterilization or disinfection medium to contact the surface of the device.
- NOTE:** Disinfection solutions do not effectively penetrate residual organic matter.
3. Prepare a mild, neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or deionized (DI) water. Immerse disassembled device and soak for 20 minutes. Use a soft brush to remove all residual organic matter, blood and irrigation solution from the device.
 4. Prepare a fresh neutral-pH enzymatic detergent (i.e., Enzol using 1 ounce per gallon) and warm (30-43°C) tap or DI water. Immerse disassembled device and perform ultrasonic cleaning for a minimum of 20 minutes.
 5. Rinse thoroughly with clean running tap or DI water.
 6. Dry all exterior surfaces with a soft cloth.
 7. Dry all openings with an air pistol (compressed air).
 8. Prior to disinfection/sterilization and reuse, this product should be inspected and reassembled to verify that it functions properly. Grasping edges and all surfaces should be smooth and free of abrasions. The device must be disassembled again prior to disinfection and/or sterilization.

NOTE: Automated cleaning methods are not compatible with the Endoscopic Curved Needle Driver.

Disinfection Instructions

NOTE: CIDEX® Activated Dialdehyde Solution is a high-level disinfection solution that has been found to be compatible with the Endoscopic Curved Needle Driver.

1. Disassemble the needle driver per the instructions listed within the IFU, if the device has been reassembled between cleaning and disinfection.
2. Activate CIDEX® Activated Dialdehyde Solution per the manufacturer's instructions. Test the activated solution prior to use with CIDEX® Solution Test Strips, to assure the solution is still activated.
3. Immerse disassembled device completely, filling all lumens or other hard-to-reach areas and eliminating air pockets, in CIDEX® Activated Dialdehyde Solution for a minimum of 45 minutes at 25°C.
4. Rinse thoroughly in a room temperature Sterile Water for Injection (SWFI) water bath. Rinse by

immersing device completely for a minimum of one minute. Manually flush all lumens and other hard-to-reach areas with SWFI water.

5. Agitate device under SWFI water, bring above water level, then re-immerses.
6. Repeat steps 4 and 5 two additional times using fresh SWFI for each of the three rinses.

NOTE: Automated disinfection methods are not compatible with the Endoscopic Curved Needle Driver.

Sterilization Instructions

NOTE: Ethylene Oxide Gas Sterilization is not recommended for the Endoscopic Curved Needle Driver.

This device can be sterilized by steam autoclave. Autoclave sterilization should be performed in accordance with industry standards with a validated sterilization cycle. The device should be disassembled prior to sterilization. The parameters needed to achieve sterility are a minimum of eight minutes at 132°C prevacuum in a single-pouched configuration. Autoclave temperature should not exceed 135°C.

LIMITED WARRANTY STATEMENT

Cook Urological Incorporated ("Company") warrants to the purchaser of this product (the "Product") that at the time of manufacture, the Product was manufactured and tested in accordance with Good Manufacturing Practices specified by the United States Food and Drug Administration and was true to label. In the unlikely event of a Product failure due to defects in material or workmanship within the warranty period, the Product will be repaired or replaced, at Cook's option, free of charge.

Because the Company has no control over the conditions under which this Product is used, diagnosis of the patient, anatomical or biological differences among patients, method of use or administration, or handling of the Product after it leaves the Company's possession, the Company does not warrant a good effect or against an ill effect following the Product's use. The foregoing warranty is exclusive and in lieu of all other warranties whether written, oral, or implied. There are no warranties of merchantability or fitness for a particular purpose. No representative of the Company may change any of the foregoing, and the purchaser accepts the Product subject to all terms hereof. The Company reserves the right to change or discontinue this Product without notice.

Limits and Exclusions

The warranty does not apply in the event of misuse, negligence, or accidental damage of Product.

Warranty services will be provided only to the original purchaser at the location where the Product was originally installed. The warranty may not be transferred without the prior written consent of the Company.

Warranty Period

This limited warranty begins on the date the Product ships from the Company's facilities and ends 6 months later.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

H. C. Topel: "Endoscopic Curved Needle Driver," Videotape, Department of Obstetrics and Gynecology, Lutheran General Hospital, Park Ridge, Illinois (available from Cook Urological and/or Cook Ob/Gyn).

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.