Novy Cornual Cannulation Set

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
Figure D

NOVY CORNUAL CANNULATION SET

The Novy Cornual Cannulation Set is intended for use through the operating channel of a hysteroscope or other uterine access device, for hysteroscopic or fluoroscopic selective catheterization and cannulation of the proximal fallopian tube(s), followed by the introduction of chromotubation solution or contrast medium, in the evaluation of tubal patency. Intended for one-time use.

The wire guide in the Novy Cornual Cannulation Set is intended only to facilitate placement of the inner catheter. It is not intended for tubal recanalization and should not be advanced beyond the tubal isthmus.

WARNING: This device should not be used in the presence of hemorrhage, active pelvic infection, profuse bleeding, or pregnancy.

NOTE: The Novy Cornual Cannulation Set is not intended for complete catheterization of the fallopian tube or tubal recanalization.

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

INDICATION STATEMENT

Selective catheterization/cannulation of the fallopian tubes for injection of dye or contrast medium in order to evaluate proximal tubal occlusion and/or patency under fluoroscopy or hysteroscopy/laparoscopy.

WARNING STATEMENTS

- **This catheter system is only designed to access the fallopian tube from the uterine cavity. It should only be advanced as far as necessary to achieve a meaningful diagnostic study. Possible adverse effects of more distal advancement have not yet been determined. The catheter system is not intended for tubal recanalization, and there is no data available, to date, showing any clinical benefit when used this way.**
- **This catheter system is not designed nor intended for embryo transfer or other assisted reproduction procedures.**
- **This device should only be used by an experienced clinician.**

CONTRAINDICATIONS

- Pregnancy or suspected pregnancy
- Active or recent pelvic infection
- Unexplained or severe uterine bleeding
- Gynecologic malignancy (suspected or known)
- Known allergy to dye or contrast medium

POSSIBLE ADVERSE EFFECTS

The following adverse effects have been reported or are of concern with the use of this type of device. Patients must be advised of the following risks.

1. **Ectopic Pregnancy.** Patients who have previously blocked fallopian tubes may be at increased risk for ectopic (tubal) pregnancy after the tubes are opened. The magnitude of the increased risk is unknown.
2. **Damage to Normal Tubes.** It is unknown whether this procedure may damage normal fallopian tubes or put the patient who has normal tubes at any increased risk for subsequent ectopic (tubal) pregnancy.
3. **Tubal Dissection.** It is possible that this catheter may dissect into the wall of the tube, under the epithelial lining and cause blockage or narrowing of the lumen of the fallopian tube. This may result in sterility, infertility, or increased risk of ectopic (tubal) pregnancy.
4. **Pain or Discomfort.** Some patients may experience cramping or other abdominal pain.

5. Extravasation. Dye or contrast medium that exceeds normal anatomic constraints of the tube may be indicative of extravasation and consequent tubal damage.
6. Tubal Perforation, Bleeding, Infection. There is a small but finite risk of tubal perforation, bleeding, and infection associated with the use of this type of device. The magnitude of this risk has not been quantified.

SET CONTAINS

- Introducing catheter with side-arm fitting and Tuohy-Borst adapter
- Stiffening obturator (*not included in all sets*)
- Inner catheter with 1 cm markings and Tuohy-Borst adapter
- Wire guide with removable silicone safety cap
- Luer lock cap

SUGGESTED INSTRUCTIONS FOR USING NOVY CORNUAL CANNULATION SET WITH HYSTEROSCOPY

NOTE: The Novy Cornual Cannulation Set should be loaded into the hysteroscope before beginning the examination, assuring the introducing catheter may be moved freely within the hysteroscope.

1. Advance the introducing catheter, with the obturator (if included) in place, into the operating channel of the hysteroscope.
2. Remove the stiffening obturator from the introducing catheter and occlude the straight arm of the fitting with the Luer lock cap provided. This port is used for irrigation or suction and injection of chromotubation solution or contrast medium into the uterine cavity.
3. Loosen the Tuohy-Borst adapter on the inner catheter and introduce the wire guide into it, positioning the distal tip of the wire guide evenly with the tip of the catheter. Tighten the Tuohy-Borst adapter to fix the position of the wire guide.
4. Loosen the Tuohy-Borst adapter on the side-arm of the introducing catheter fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until both catheter tips are even. For future reference, note the position of the proximal distance markings on the inner catheter with respect to the Tuohy-Borst adapter on the side-arm fitting. The first marking should be flush with the Tuohy-Borst adapter.
5. Remove the inner catheter and wire guide from the introducing catheter, and tighten the Tuohy-Borst adapter on the side-arm fitting to seal it and prevent leakage of fluid or gas.
6. Position the introducing catheter in the operating channel of the hysteroscope so it does not extend beyond the telescope lens.
7. Introduce the hysteroscope sheath, with its blunt obturator in place, through the dilated cervix and into the uterine cavity.
8. Remove the obturator from the hysteroscope sheath.
9. Introduce the hysteroscope, with the introducing catheter in place, into the hysteroscope sheath.
10. Under direct vision, guide the hysteroscope toward the tubal ostium.
11. After the tubal ostium is visualized, direct the introducing catheter toward the tubal orifice and wedge the tip of the catheter in the tubal ostium.
12. Remove the Luer lock cap from the straight arm of the side-arm fitting and attach a 3 ml syringe filled with chromotubation solution or dilute contrast medium (25% solution).
13. Inject the chromotubation solution or contrast medium into the tubal ostium while observing laparoscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency has been confirmed and the procedure is completed.

NOTE: If spillage of chromotubation solution or contrast medium from the fimbrial portion of the tube into the peritoneal cavity is not observed, it will be necessary to introduce the inner catheter into the proximal fallopian tube to confirm tubal blockage or rule-out a false positive result due to temporary tubal spasm. Please proceed with step 14.

14. Detach the 3 ml syringe and occlude the straight arm of the side-arm fitting with the Luer lock cap.
15. Loosen the Tuohy-Borst adapter on the side-arm fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until the first proximal marking on the inner catheter is flush with the Tuohy-Borst adapter. The distal tip of the inner catheter should now be even with that of the introducing catheter.
16. Taking care to allow movement of the inner catheter, tighten the Tuohy-Borst adapter on the side-arm fitting to form a gas tight seal around the inner catheter.
17. Cannulation of the intramural portion of the fallopian tube with the wire guide and inner catheter is begun by first advancing the wire guide (Figure A). The flexibility of the wire guide tip will increase as the length exiting the catheter tip increases.
18. Advance the wire guide beyond the utero-tubal junction and into the tubal isthmus unless significant resistance is met (Figure B). Advancement of the wire guide through the utero-tubal junction is observed directly by laparoscopy and video display of the proximal oviduct or by characteristic arching of the tubal isthmus. Advancement may be facilitated by displacing or straightening the fallopian tube with an instrument under laparoscopic guidance.

NOTE: If resistance is met and it is not possible to easily advance the wire guide into the tubal isthmus, exchange the wire guide included in the set for a .015 inch wire guide (order numbers: J-SMGC-150601 or J-SMGC-150901) and try again.

19. Coaxially advance the inner catheter over the positioned wire guide for a short distance or until resistance is felt (Figures C and D). The distance the catheter passes into the tubal isthmus may be determined by observing the distance markings on either the proximal or distal portion of the catheter.

WARNING: If significant resistance is felt, do not attempt to advance the catheter.

WARNING: Advancement of the catheter or wire guide beyond the tubal isthmus should not be attempted.

20. Remove the wire guide and attach a 3 ml syringe filled with chromotubation solution or dilute contrast medium (25% solution).
21. Inject the chromotubation solution or contrast medium into the tubal ostium while observing laparoscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency is confirmed.

SUGGESTED INSTRUCTIONS FOR USING NOVY CORNUAL CANNULATION SET WITH FLUOROSCOPY

1. Remove the stiffening obturator (if included) from the introducing catheter and occlude the straight arm of the fitting with the Luer lock cap provided. This port is used for irrigation or suction and injection of contrast medium into the uterine cavity.
2. Loosen the Tuohy-Borst adapter on the inner catheter and introduce the wire guide into it, positioning the distal tip of the wire guide evenly with the tip of the catheter. Tighten the Tuohy-Borst adapter to fix the position of the wire guide.
3. Loosen the Tuohy-Borst adapter on the side-arm of the introducing catheter fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until both catheter tips are even. For future reference, note the position of the proximal distance markings on the inner catheter with respect to the Tuohy-Borst adapter on the side-arm fitting. The first marking should be flush with the Tuohy-Borst adapter.
4. REMOVE the inner catheter and wire guide from the introducing catheter, and tighten the Tuohy-Borst adapter on the side-arm fitting to seal it and prevent leakage.

5. Introduce the appropriate uterine access device into the cervix and affix it in position either on the cervix or within the uterine cavity.
6. Insert the introducing catheter through the uterine access device and advance the catheter until its distal tip is positioned in the lower portion of the uterine cavity.
7. Remove the Luer lock cap from the straight arm of the side-arm fitting and attach a syringe filled with half-strength contrast medium. With fluoroscopic visualization, inject contrast medium through the introducing catheter to aid visualization of the uterine cavity and tubal ostium.
8. Detach the syringe and introduce a .035 inch Safe-T-J® wire guide (order number: J-SCF-350903) through the introducing catheter and into the uterine cavity. With the wire guide leading and under fluoroscopic guidance, advance the wire guide and catheter together toward the tubal ostium, and wedge the tip of the catheter in the tubal ostium.
9. Remove the wire guide from the introducing catheter and attach a syringe filled with full-strength contrast medium. Inject 2 to 10 ml contrast medium into the tubal ostium while observing fluoroscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency has been confirmed and the procedure is completed.

NOTE: If spillage of contrast medium from the fimbrial portion of the tube into the peritoneal cavity is not observed, it will be necessary to introduce the inner catheter into the proximal fallopian tube to confirm tubal blockage or rule-out a false positive result due to temporary tubal spasm. Please proceed with step 10.

10. Detach the syringe and occlude the straight arm of the side-arm fitting with the Luer lock cap.
11. Loosen the Tuohy-Borst adapter on the side-arm fitting and advance the inner catheter, with the wire guide in place, through the introducing catheter until the first proximal marking on the inner catheter is flush with the Tuohy-Borst adapter. The distal tip of the inner catheter should now be even with that of the introducing catheter.
12. Taking care to allow movement of the inner catheter, tighten the Tuohy-Borst adapter on the side-arm fitting to form a fluid tight seal around the inner catheter.
13. Cannulation of the intramural portion of the fallopian tube with the wire guide and inner catheter is begun by first advancing the wire guide (Figure A). The flexibility of the wire guide tip will increase as the length exiting the catheter tip increases.
14. Advance the wire guide beyond the utero-tubal junction and into the tubal isthmus unless significant resistance is met (Figure B). Advancement of the wire guide through the utero-tubal junction may be observed fluoroscopically.

NOTE: If resistance is met and it is not possible to easily advance the wire guide into the tubal isthmus, exchange the wire guide included in the set for a .015 inch wire guide (order numbers: J-SMGC-150601 or J-SMGC-150901) and try again.

15. Coaxially advance the inner catheter over the positioned wire guide for a short distance or until resistance is felt (Figures C and D). The distance the catheter passes into the tubal isthmus may be determined by observing the distance markings on the proximal portion of the catheter.

WARNING: If significant resistance is felt, do not attempt to advance the catheter.

WARNING: Advancement of the catheter or wire guide beyond the tubal isthmus should not be attempted.

16. Remove the wire guide and attach a syringe filled with full-strength contrast medium.
17. Inject the contrast medium into the tubal ostium while observing fluoroscopically for spillage from the fimbrial portion of the tube. If such spillage is observed, tubal patency is confirmed.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

REFERENCE

M. J. Novy, A. S. Thurmond, P. Patton, B. T. Uchida, J. Rösch: "Diagnosis of Cornual Obstruction by Transcervical Fallopian Tube Cannulation," *Fertility and Sterility*, 50 (1988), 434-440.5.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.