EN
3

Novapass™ Suture Passers

Instructions for Use

DA
4

Novapass™ suturpassere

Brugsanvisning

DE
5

Novapass™ Fadentransporter

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EL
6

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




5

NOVAPASS™ SUTURE PASSERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner)

DEVICE DESCRIPTION

Novapass Suture Passers are hand-held suture grasping instruments designed to pass suture through soft tissue. Novapass Suture Passers are available in three configurations, as shown in **Figs. 1-3** and detailed in the table below.

	Gage	Length	For Use in:
30° Reverdin Curve (Fig. 1) Part number: C-PAS-30-11.4	12 gage	11.4 cm	Open abdominal procedures
60° Reverdin Curve (Fig. 2) Part Number: C-PAS-60-11.8	12 gage	11.8 cm	Open abdominal procedures
Straight Suture Passer (Fig. 3) Part Number: C-PAS-180-16.0	12 gage	16.0 cm	Laparoscopic or open abdominal procedures

INTENDED USE

Novapass Suture Passers are intended to pass suture through soft tissue for open and laparoscopic abdominal procedures.

CONTRAINDICATIONS

None known

WARNINGS

Exercise extreme caution during use of this device to avoid unintentional puncture/perforation of vessel, nerve, bowel, or any other internal organs or vital anatomical structures.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in open and laparoscopic surgical techniques. Standard techniques should be employed.
- Do not resterilize this device. Discard all open and unused portions of the suture passer.
- **Do not attempt to repair a bent or broken suture passer.** Discard the device if mishandling has caused possible damage or contamination. Use of a damaged product may result in failure of the device.
- Do not use the suture passer in a procedure where the position of the needle tip cannot be clearly identified.
- Novapass Suture Passers have needle-point tips that can cause injury to internal tissue and/or any personnel who may come in contact with the device. The needle-point tip should be kept protected with the protective cap whenever the device is not in use.
- Novapass Suture Passers are recommended for use only with suture sizes of 2, 1, 0, 2-0, and 3-0. Use of the device with smaller or larger size sutures may damage the suture or prevent the device from functioning properly.
- The needle-point tip of the device must be in the fully retracted/closed position prior to passing in or out of tissue. Maintain pressure on the trigger or keep the trigger in the "LOCKED" (**Fig. 4**) position to ensure full retraction of the needle-point tip during passage in and out of tissue.

INSTRUCTIONS FOR USE

1. Remove the protective cap from the needle-point tip of the device. **NOTE:** Take care to avoid user injury caused by the needle-point tip. Retain the protective cap to protect the needle-point tip when the device is not in use.
2. Push the trigger forward into the "OPEN" (**Fig. 5**) position to expose the suture hook.
3. Place the operative suture into the hook.
4. Pull back on the trigger to secure the suture within the hook. The trigger can also be locked (**Fig. 4**) to prevent the suture from inadvertently slipping out of the device.
5. Prior to passing the device in or out of tissue, ensure the needle tip is fully closed/retracted by pulling back on the trigger away from the "OPEN" (**Fig. 5**) position.
6. Pass the device through the desired tissue. Take care to avoid unintentional puncture/perforation of vessel, nerve, bowel, and any other internal organs or vital anatomical structures.
7. Release or unlock the trigger and push forward toward the "OPEN" (**Fig. 5**) position to retrieve the suture.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.