**EN**  
**3**

# **Rutner Percutaneous Suprapubic Balloon Catheter Set**

Instructions for Use

**CS**  
**4**

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## RUTNER PERCUTANEOUS SUPRAPUBIC BALLOON CATHETER SET

**CAUTION:** U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

### DEVICE DESCRIPTION

- Polyurethane balloon catheter, 22 cm long
- Stainless steel needle obturator
- Clear polyvinylchloride connecting tube with one-way stopcock and drainage bag connector
- Syringe

### INTENDED USE

Used to provide bladder drainage by percutaneous placement of a balloon catheter.

### CONTRAINDICATIONS

None known

### WARNINGS

Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.

### PRECAUTIONS

- Catheter must not remain indwelling more than four weeks.
- This product contains **natural rubber latex** which may cause allergic reactions in some individuals.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### INSTRUCTIONS FOR USE

1. Remove the protective sleeve from the balloon and discard. Place the needle obturator inside the balloon catheter and secure its position with the Luer lock.
2. Fill the bladder with sterile saline solution if not already distended with urine.
3. Make a 3 to 4 mm incision in the skin above the symphysis pubis, over or just below the palpable bladder dome.

**CAUTION:** Using excessive pressure to inflate the latex balloon on this device can cause the balloon to rupture. Inflation of the balloon 2-3 times with an appropriate inflation medium prior to use will allow easier inflation of the balloon.

4. Insert the balloon catheter, using one hand to push the catheter from behind the needle hub (using palm of hand), and the other at skin level to guide the needle. Insert the balloon catheter vertically into the bladder on patients without previous pelvic surgery or at a 30° angle toward the symphysis pubis on patients with previous pelvic surgery. Advance the catheter no more than 4-5 cm beyond the first sign of bladder fluid coming out of the needle obturator. The balloon catheter midpoint (etch mark on the needle) is an approximation.
5. If bladder pressure is low and urine fails to flow spontaneously through the needle obturator, aspirate urine to ascertain that the balloon catheter is within the bladder.

**WARNING:** Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide or any other gas.

6. Fill the syringe with sterile water or saline solution. Eliminate air from the syringe.
7. Inject appropriate amount of sterile water or saline solution into the balloon catheter. **CAUTION:** Do not overinflate the balloon. Refer to product label or the inflation check valve on the balloon device for appropriate balloon volume.
8. Release and remove the needle obturator.
9. Attach the connecting tube and one-way stopcock to the balloon catheter and secure to a drainage bag.
10. Allow the bladder to empty.

11. Withdraw the balloon catheter until the balloon meets the resistance of the bladder dome (**Fig. 1**).
  12. Apply tincture of benzoin to the balloon catheter shaft and the abdominal skin.
  13. Secure the catheter to the skin with a strip of 2 inch wide tape.
  14. Wrap the tape completely around the balloon catheter shaft, forming a 1 cm adhesive land to allow for catheter movement; lay the tape flat on the abdomen.
- NOTE:** Skin sutures are unnecessary for retention of this catheter.
15. Secure the connecting tube to the abdomen in a gentle curve with tape (**Fig. 2**).
  16. Check to ensure that the one-way stopcock handle is in the direction of flow for urine drainage through the catheter. Turn handle 90° to the direction of drainage to obstruct the flow when appropriate for trials of voiding.

**NOTE:** It is recommended the balloon's fluid volume be checked every 3-5 days and then reinflated with recommended volume.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.